TO GUIDE AND FOCUS LIGHT AND TO SELECT A WAVELENGTH, THE SELECT SOURCE USES MIRRORS, NETWORKS, AND LENSES, AS IN THE VISIBLE SPECTRUM. BUT THE MATERIALS AND OPTICAL CONDITIONS VARY ACCORDING TO THE WAVELENGTHS USED. THAT IS WHY SOLEIL BEAMLINES ARE EQUIPPED WITH A WIDE VARIETY OF OPTICS.

THICK STEEL ENCLOSURES

Because the beam propagates in a vacuum, the slightest deformation of the

mechanisms can result in its misalignment. To ensure adequate stability of the line,

the optical devices are placed in thick stainless steel enclosures.

OPTICAL DEVICES TO 'PREPARE' THE BEAM

In each beamline, the optical instruments give the beam the characteristics that are required for the experiment. They direct it, select a wavelength, and concentrate the greatest possible number of photons on the sample.

EXAMPLE OF A BEAMLINE IN THE X-RAY DOMAIN: CRYSTAL MONOCHROMATOR s continuously adjusted as the angular position of the crystals is varied (Bragg diffraction). To cover a more complete Two perpendicular slits are inenergy domain, different crystals are terposed on the optical path to **Optics hutch** The filter absorbs low-energy photons, mounted on water-cooled supports delineate the light beam (only particularly infra-red and visible. It is a 1 and 2 : silicon 111, beryllium, etc. one is shown here). These sheet of aluminium, carbon, or copper, Each support can contain up to five slits can be in copper, steel or or even a diamond, which eliminates crystals. They are selected by simply tungsten, but any x-ray absorapplying a translation to the device. some of the power of the beam. This bent material is suitable Workstation prevents the overheating of the optics positioned downstream. Effect exploited: selective absorption The monochromator selects a type of monochromatic light from The two crystal supports the initial 'white' beam. It is an assembly of two silicon crystals. A computer controls their angular position to vary the wavelength and select the one that will make the sample react. Effect exploited: Bragg diffraction. To focus light on the sample, there are two consecutive mirrors whose curvatures are perpendicular, in a Kirkpatrick-Baez (KB) configuration. The shape of the curvatures is computer-controlled in order to adjust the focal Slit Effect exploited: total reflection. The sample reacts when it is lit by the selected monochromatic light. The experiment consists of recording Filter the response using a detector suitable for the type of signal emitted (absorption spectroscopy, fluorescent emission, KB mirror enclosure diffraction, etc.). Monochromator The pair of perpendicular KB mirrors KIRKPATRICK-BAEZ (KB) MIRRORS Mirrors The Kirkpatrick-Baez (KB) assembly consists of two perpendicular mirrors 3 & 4 whose longitudinal profiles are adjusted by computer-controlled 'benders'. By applying bending 7.5 m 3.5 m 4.5 m 1 m forces to each end, the benders can give an elliptical shape to flat Detector — Sample or toric mirrors. ph/s/mm² (1.24 10¹⁶